

# Adjustable instrument table

## AIT-16

- Motorized height adjustment
- Smooth castor movement with locking mechanism
- Load capacity 50kg

The TOPCON AIT-16 is specifically designed for mounting tabletops of TOPCON ophthalmic and medical instruments. With a maximum load capacity of 50kg, even a complete retinal camera can be accommodated.

- The up-down switch is conveniently located on the column of the table.
- The compact style allows placing instruments on a small footprint.
- Two auxiliary power outlets are provided.
- A tabletop is available as an option.

### Specifications

Maximum table height (w/Tabletop)	880mm
Minimum table height (w/Tabletop)	660mm
Vertical adjustment range	220mm
Maximum load	50kg
Overall dimensions (body only)	510(W) x 450(D) x 640~860(H)mm
Weight (body only)	23kg
Power supply	AC 230 Volts
Power consumption	270 VA
Frequency	50/60 Hz



TOPCON CORPORATION

TOPCON CORPORATION

\*Subject to change in design and/or specifications without advanced notice.

#### IMPORTANT

In order to obtain the best results with this instrument, please be sure to review all user instructions prior to operation.